

CLUSTER ANALYSIS


Visual instructions by Laura Martin, February 2022 1st 2022

Methods in
Molecular Biology 1480

Springer Protocols

STORM microscopy and cluster
analysis for PcG studies

Laura Martin, Álvaro Castells-García, Maria Pia Cosma,
and Maria Victoria Neguembor




CLUSTER ANALYSIS

Original script from Carlo Manzo

Modified in the years by
J Otterstorm, J Borbely, A Castells, & me

20 functions



CLUSTER ANALYSIS

Hierarchical code organization

>CLUSTER_ANALYSIS_main

>defineClustersStruct
>FindClusterStruct



Set analysis **parameters**

>FindClusters



Identifies Islands and clusters
Calculates clusters' features:

- n° localizations
- Areas
- NND (of In-Island only)

>saveClusterMetricData
>extractClusterStats
>DistanceDualColor



Calculates the **global NND**
between ALL clusters

Saves output files
(.xlsx, .mat, .bin,
.dcc, .fig, .png)

CLUSTER ANALYSIS

Input Parameters for Find...

image_width [pixels]
256

image_height [pixels]
256

use_drift_corrected_xy (1/0 true/false)
1

use_channels [all channels: -1]
-1

original_pixel_size [nm]
160

analysis_pixel_size [nm]
10

sum_roi_size [pixels]
5

sum_threshold
5

factor [use: analysis_pixel / factor]
5

localization_precision [nm]
9

minimum_molecules_per_cluster
5

ignoreNumPeakThreshold
20000

ignoreNumIslandThreshold
15000

drawUsingZRange (1/0 true/false)
0

use_iterative_segmentation (1/0 true/false)
1

max_segmentation_area
20000

show_density_map (1/0 true/false)
0

show_mask (1/0 true/false)
0




OK Cancel

Dependent on
Microscope/Camera
settings

CRUCIAL FOR CLUSTER
ANALYSIS


Optionals
(Particular conditions, plots...)

CLUSTER ANALYSIS



CLUSTER ANALYSIS

10 nm




original_pixel_size [nm]

160

analysis_pixel_size [nm]

10

CLUSTER ANALYSIS



original_pixel_size [nm]

160

analysis_pixel_size [nm]

10

sum_roi_size [pixels]


5

sum_threshold

5

CLUSTER ANALYSIS


10 nm



ISLANDs identification

CLUSTER ANALYSIS


10 nm



It picks one Island


CLUSTER ANALYSIS

10 nm
┌───┐




It picks one Island

CLUSTER ANALYSIS



It picks one Island

It resizes it



analysis_pixel_size [nm]


10

factor [use: analysis_pixel / factor]

5

10 nm : **5** = 2 nm


CLUSTER ANALYSIS




It picks one Island

It resizes it

For each
localization:




CLUSTER ANALYSIS



$$\longrightarrow p(x) = \frac{1}{\sigma\sqrt{2\pi}} e^{-\frac{1}{2}\left(\frac{x-\mu}{\sigma}\right)^2}$$


It picks one Island

It resizes it





For each
localization:
It associates
a Gaussian curve


CLUSTER ANALYSIS




$$\longrightarrow p(x) = \frac{1}{\sigma\sqrt{2\pi}} e^{-\frac{1}{2}\left(\frac{x-\mu}{\sigma}\right)^2}$$



x, y localization coordinates





$$\sigma = \frac{\text{loc precision (9nm)}}{\text{new analysis px (2nm)}}$$



$$\text{roi size} = \text{numSigma} * \sigma$$

CLUSTER ANALYSIS



New matrix of the
summation of
Gaussian values
belonging to
localizations


It picks one Island

It resizes it

For each
localization:
It associates
a Gaussian curve

It sums the
Gaussians to
generate a
density map

CLUSTER ANALYSIS



In each Island,
from the density matrix
the Matlab function **regionprops**
'magically' identifies centroids
of clusters

Then, the function **lik_sig** retrieves the Number of localizations belonging to each centroid, and the coordinates of the centroids in nm.

If the Number of localizations
belonging to each centroid is
 \geq `minimum_molecules_per_cluster`

a cluster is defined.

`minimum_molecules_per_cluster`

5

CLUSTER ANALYSIS



N° of localizations per cluster are:
the sum of localizations belonging to the cluster

For each cluster, the algorithm calculates the
coordinates **X, Y of the CENTROID** as:
the mean of localizations coordinates (x, y)

CLUSTER ANALYSIS

Area of the cluster is approximated to a circle, of which

$$r = \frac{sd\ X + sd\ Y}{2} \quad A = \pi r^2$$

*To convert px to nm, simply *160 (original pixel size)*

Inside each Island, distance among all cluster centroids (X, Y) is calculated, and the smallest distance between each pair of clusters is given in nm as **NND**

$$distance = \sqrt{(X_b - X_a)^2 + (Y_b - Y_a)^2}$$

Of course, for Single clusters, NND will be 'Infinite' since for definition they have no close neighbours

CLUSTER ANALYSIS

>**DistanceDualColor**



To calculate NND
among ALL clusters

NND is calculated using Centroids coordinates (X, Y)
from EVERY cluster (Inland and Single)

This script was originally made to calculate distances
among different proteins (dual color).

For single color, the list of centroids is “doubled” and
distances are calculated between the two identical
lists. (distances = 0 are discarded.)